

# How to Use Fentanyl Test Strips<sup>1</sup>

## 1 Prepare

- Take a 5-10 mg scoop of the substance
- Add substance to about 1 teaspoon of clean water, or enough to fill the cap of a water bottle
- Mix the substance and water in a dry container (a cooker, bottle cap, or baggie the substance came in)




<sup>2</sup> 10 mg = enough to cover the hair of Abraham Lincoln's head on a penny

## 2 Test substance


- Hold the **solid blue** end of the test strip
- Dip the test strip into the water up to the red line for 15 seconds then take it out
- Wait for 1-2 minutes to check your results



## 3 Read results



**ONE red line:**  
**Positive for fentanyl**  
This means fentanyl is present



**TWO red lines**  
**Negative for fentanyl**  
This means fentanyl is not present

The second line can be faint,  
this still means negative



Scan here to visit our website for more resources and information



**MATTERS**  
Medication for Addiction Treatment & Electronic Referrals